


Product Datasheet

TIGIT antibody (orb1326412)

Description	Mouse monoclonal antibody to TIGIT
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Full length human recombinant protein of human TIGIT (NP_776160) produced in HEK293T cell.
Target	TIGIT
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	26.1 kDa
Uniprot ID	Q495A1
NCBI	NM_173799
Dilution Range	FLOW 1:100
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of living 293T ...