


Product Datasheet

Complement C4B (C4B) antibody (orb1326428)

Description	Mouse monoclonal antibody to Complement C4B (C4B)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1100-1350 of human C4B (NP_001002029) produced in E.coli.
Target	C4B
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P0C0L4
NCBI	NM_001002029
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...