


Product Datasheet

PD-L1 (CD274) antibody (orb1326443)

Description	Mouse monoclonal antibody to PD-L1 (CD274)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human CD274 (NP_054862) produced in HEK293T cell.
Target	CD274
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	33.28 kDa
Uniprot ID	Q9NZQ7
NCBI	NM_014143
Dilution Range	WB 1:500, FLOW 1:50~100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...