


Product Datasheet

Cytoglobin (CYGB) antibody (orb1326465)

Description	Mouse monoclonal antibody to Cytoglobin (CYGB)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human CYGB (NP_599030) produced in HEK293T cell.
Target	CYGB
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	21.2 kDa
Uniprot ID	Q8WWM9
NCBI	NM_134268
Dilution Range	WB 1:250~500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...