

Product Datasheet

C20orf79 (SCP2D1) antibody (orb1326480)

Description Mouse monoclonal antibody to C20orf79 (SCP2D1)

Species/Host Mouse

Reactivity Human

Conjugation Unconjugated

Tested Applications WB

Immunogen Full length human recombinant protein of human C20orf79 (NP_848578) produced in E.coli.

Target SCP2D1

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG1

Clonality Monoclonal


MW 17.5 kDa

Uniprot ID [Q9UJQ7](#)

NCBI [NM_178483](#)

Dilution Range WB 1:500

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...