


Product Datasheet

ATF 4 (ATF4) antibody (orb1326498)

Description	Mouse monoclonal antibody to ATF 4 (ATF4)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human Atf4 (NP_001666) produced in E.coli.
Target	ATF4
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	38.4 kDa
Uniprot ID	P18848
NCBI	NM_001675
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...