


Product Datasheet

PSMA1 antibody (orb1326514)

Description	Mouse monoclonal antibody to PSMA1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human PSMA1 (NP_002777) produced in E.coli.
Target	PSMA1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	29.4 kDa
Uniprot ID	P25786
NCBI	NM_002786
Dilution Range	WB 1:250
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts (35ug)...