


Product Datasheet

EFHD1 antibody (orb1326520)

Description	Mouse monoclonal antibody to EFHD1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human EFHD1 (NP_079478) produced in HEK293T cell.
Target	EFHD1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	26.7 kDa
Uniprot ID	Q9BUP0
NCBI	NM_025202
Dilution Range	WB 1:250
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts (35ug)...