


Product Datasheet

CAD antibody (orb1326531)

Description	Mouse monoclonal antibody to CAD
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1918-2225 of human CAD (NP_004332) produced in E.coli.
Target	CAD
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P27708
NCBI	NM_004341
Dilution Range	WB 1:250~500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...