

Product Datasheet

IL28 Receptor alpha (IFNLR1) antibody (orb1326586)

Description Mouse monoclonal antibody to IL28 Receptor alpha (IFNLR1)

Species/Host Mouse

Reactivity Human

Conjugation Unconjugated

Tested Applications IHC, WB

Immunogen Human recombinant protein fragment corresponding to amino acids 21-228 of human IFNLR1 (NP_734464) produced in E.coli.

Target IFNLR1

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG1


Clonality Monoclonal

Uniprot ID [Q8IU57](#)

NCBI [NM_170743](#)

Dilution Range WB 1:500, IHC 1:500

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...