


Product Datasheet

TIM 3 (HAVCR2) antibody (orb1326604)

Description	Mouse monoclonal antibody to TIM 3 (HAVCR2)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human HAVCR2 (NP_116171) produced in HEK293T cell.
Target	HAVCR2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG3
Clonality	Monoclonal
MW	33.2 kDa
Uniprot ID	Q8TDQ0
NCBI	NM_032782
Dilution Range	WB 1:500, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...