


Product Datasheet

MRP2 (ABCC2) antibody (orb1326643)

Description	Mouse monoclonal antibody to MRP2 (ABCC2)
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1256-1545 of human ABCC2 (NP_000383) produced in E.coli.
Target	ABCC2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q92887
NCBI	NM_000392
Dilution Range	WB 1:500, IHC 1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...