

Product Datasheet

LDL Receptor (LDLR) antibody (orb1326651)

Description Mouse monoclonal antibody to LDL Receptor (LDLR)

Species/Host Mouse

Reactivity Human

Conjugation Unconjugated

Tested Applications FC, WB

Immunogen Full length human recombinant protein of human LDLR (NP_000518) produced in HEK293T cell.

Target LDLR

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG1

Clonality Monoclonal


MW 95.38 kDa

Uniprot ID [P01130](#)

NCBI [NM_000527](#)

Dilution Range WB 1:500~2000, FLOW 1:100

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...