


Product Datasheet

DIP13B (APPL2) antibody (orb1326690)

Description	Mouse monoclonal antibody to DIP13B (APPL2)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human APPL2 (NP_060641) produced in HEK293T cell.
Target	APPL2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	74.3 kDa
Uniprot ID	Q8NEU8
NCBI	NM_018171
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...