


Product Datasheet

FOXN2 antibody (orb1326708)

Description	Mouse monoclonal antibody to FOXN2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-431 of human FOXN2 (NP_002149) produced in E.coli.
Target	FOXN2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	47 kDa
Uniprot ID	P32314
NCBI	NM_002158
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...