


Product Datasheet

KAT5 antibody (orb1326719)

Description	Mouse monoclonal antibody to KAT5
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 227-504 of human KAT5 (NP_006379) produced in E.coli.
Target	KAT5
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	58.4 kDa
Uniprot ID	Q92993
NCBI	NM_006388
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts (35ug)...