


Product Datasheet

GI24 (C10orf54) antibody (orb1326753)

Description	Mouse monoclonal antibody to GI24 (C10orf54)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human C10ORF54 (NP_071436) produced in HEK293T cell.
Target	VSIR
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	33.7 kDa
Uniprot ID	Q9H7M9
NCBI	NM_022153
Dilution Range	WB 1:500, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...