


Product Datasheet

IL15RA antibody (orb1326767)

Description	Mouse monoclonal antibody to IL15RA
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 31-267 of human IL15RA (NP_001230468) produced in E.coli.
Target	IL15RA
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	28.22 kDa
Uniprot ID	Q13261
NCBI	NM_001243539
Dilution Range	WB 1:500, FLOW 1:25
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...