


Product Datasheet

LYVE1 antibody (orb1326771)

Description	Mouse monoclonal antibody to LYVE1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Human recombinant protein fragment corresponding to amino acids 20-238 of human LYVE1 (NP_006682) produced in HEK293T cell.
Target	LYVE1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	35 kDa
Uniprot ID	Q9Y5Y7
NCBI	NM_006691
Dilution Range	FLOW 1:100
Expiration Date	12 months from date of receipt.



Flow
cytometric
Flow
cytometric
HEK293T ...