


Product Datasheet

Collagen III (COL3A1) antibody (orb1326809)

Description	Mouse monoclonal antibody to Collagen III (COL3A1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 981-1205 of human COL3A1 (NP_000081) produced in E.coli.
Target	COL3A1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	138.56 kDa
Uniprot ID	P02461
NCBI	NM_000090
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...