


Product Datasheet

CD11d (ITGAD) antibody (orb1326834)

Description	Mouse monoclonal antibody to CD11d (ITGAD)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 18-418 of human ITGAD (NP_005344) produced in E.coli.
Target	ITGAD
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	126.76 kDa
Uniprot ID	Q13349
NCBI	NM_005353
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...