


Product Datasheet

CREM antibody (orb1326845)

Description	Mouse monoclonal antibody to CREM
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human CREM (NP_874388) produced in E.coli.
Target	CREM
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	13.3 kDa
Uniprot ID	Q03060
NCBI	NM_182719
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...