


Product Datasheet

Cornulin (CRNN) antibody (orb1326971)

Description	Mouse monoclonal antibody to Cornulin (CRNN)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human CRNN (NP_057274) produced in E.coli.
Target	CRNN
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	53.4 kDa
Uniprot ID	Q9UBG3
NCBI	NM_016190
Dilution Range	WB 1:2000, IHC 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...