


Product Datasheet

Factor X (F10) antibody (orb1327031)

Description	Mouse monoclonal antibody to Factor X (F10)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 41-488 of human F10 (NP_000495) produced in E.coli.
Target	F10
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	54.54 kDa
Uniprot ID	P00742
NCBI	NM_000504
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...