


Product Datasheet

PDPN antibody (orb1327061)

Description	Mouse monoclonal antibody to PDPN
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 23-131 of human PDPN (NP_006465) produced in E.coli.
Target	PDPN
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	24.7 kDa
Uniprot ID	Q86YL7
NCBI	NM_006474
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...