


Product Datasheet

ZFP57 antibody (orb1327117)

Description	Mouse monoclonal antibody to ZFP57
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human ZFP57 (NP_001103279) produced in E.coli.
Target	ZFP57
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q9NU63
NCBI	NM_001109809
Dilution Range	WB 1:500~2000, IHC 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...