


Product Datasheet

Repulsive Guidance Molecule A (RGMA) antibody (orb1327128)

Description	Mouse monoclonal antibody to Repulsive Guidance Molecule A (RGMA)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Human recombinant protein fragment corresponding to amino acids 169-424 of human RGMA (NP_064596) produced in E.coli.
Target	RGMA
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	49.2 kDa
Uniprot ID	Q96B86
NCBI	NM_020211
Dilution Range	IHC 1:250
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...