


Product Datasheet

FAM119A (METTL21A) antibody (orb1327140)

Description	Mouse monoclonal antibody to FAM119A (METTL21A)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human METTL21A (NP_001120867) produced in E.coli.
Target	METTL21A
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	24.4 kDa
Uniprot ID	Q8WXB1
NCBI	NM_001127395
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...