


Product Datasheet

SAPAP3 (DLGAP3) antibody (orb1327158)

Description	Mouse monoclonal antibody to SAPAP3 (DLGAP3)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 794-902 of human DLGAP3 (NP_001073887) produced in E.coli.
Target	DLGAP3
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	105.9 kDa
Uniprot ID	O95886
NCBI	NM_001080418
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...