


Product Datasheet

N WASP (WASL) antibody (orb1327285)

Description	Mouse monoclonal antibody to N WASP (WASL)
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human WASL (NP_003932) produced in HEK293T cell.
Target	WASL
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	54.6 kDa
Uniprot ID	O00401
NCBI	NM_003941
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...