


Product Datasheet

NAGA antibody (orb1327291)

Description	Mouse monoclonal antibody to NAGA
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human NAGA (NP_000253) produced in HEK293T cell.
Target	NAGA
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	46.56 kDa
Uniprot ID	P17050
NCBI	NM_000262
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...