


Product Datasheet

ZNF19 antibody (orb1327298)

Description	Mouse monoclonal antibody to ZNF19
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human ZNF19 (NP_008892) produced in E.coli.
Target	ZNF19
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P17023
NCBI	NM_006961
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...