


Product Datasheet

CCDC27 antibody (orb1327326)

Description	Mouse monoclonal antibody to CCDC27
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human CCDC27 (NP_689705) produced in E.coli.
Target	CCDC27
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	75.2 kDa
Uniprot ID	Q2M243
NCBI	NM_152492
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...