


Product Datasheet

AMCase (CHIA) antibody (orb1327327)

Description	Mouse monoclonal antibody to AMCase (CHIA)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human CHIA (NP_068569) produced in E.coli.
Target	CHIA
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	40 kDa
Uniprot ID	Q9BZP6
NCBI	NM_021797
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...