


Product Datasheet

SF3B14 (SF3B6) antibody (orb1327388)

Description	Mouse monoclonal antibody to SF3B14 (SF3B6)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human SF3B14 (NP_057131) produced in E.coli.
Target	SF3B6
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	14.4 kDa
Uniprot ID	Q9Y3B4
NCBI	NM_016047
Dilution Range	WB 1:2000, IHC 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...