


Product Datasheet

CIB1 antibody (orb1327389)

Description	Mouse monoclonal antibody to CIB1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 2-191 of human CIB1 (NP_006375) produced in E.coli.
Target	CIB1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	21.5 kDa
Uniprot ID	Q99828
NCBI	NM_006384
Dilution Range	WB 1:500~2000, IHC 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...