

Product Datasheet

CSPS (SULT1A3) antibody (orb1327405)

Description	Mouse monoclonal antibody to CSPS (SULT1A3)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human SULT1A3 (NP_003157) produced in E.coli.
Target	SULT1A3
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	34 kDa
Uniprot ID	P50224
NCBI	NM_003166
Dilution Range	WB 1:500~2000, IHC 1:500
Expiration Date	12 months from date of receipt.

