


Product Datasheet

ZNF202 antibody (orb1327468)

Description	Mouse monoclonal antibody to ZNF202
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-307 of human ZNF202 (NP_003446) produced in E.coli.
Target	ZNF202
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	74.5 kDa
Uniprot ID	O95125
NCBI	NM_003455
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...