


Product Datasheet

ASH2L antibody (orb1327490)

Description	Mouse monoclonal antibody to ASH2L
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 347-628 of human ASH2L (NP_004665) produced in E.coli.
Target	ASH2L
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	68.5 kDa
Uniprot ID	Q9UBL3
NCBI	NM_004674
Dilution Range	WB 1:2000, IHC 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...