


Product Datasheet

JMJD5 (KDM8) antibody (orb1327516)

Description	Mouse monoclonal antibody to JMJD5 (KDM8)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human JMJD5 (NP_079049) produced in HEK293T cell.
Target	KDM8
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	47.1 kDa
Uniprot ID	Q8N371
NCBI	NM_024773
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...