

Product Datasheet

HNRPH1 (HNRNPH1) antibody (orb1327558)

Description Mouse monoclonal antibody to HNRPH1 (HNRNPH1)

Species/Host Mouse

Reactivity Human, Mouse

Conjugation Unconjugated

Tested Applications IHC, WB

Immunogen Full length human recombinant protein of human HNRNPH1 (NP_005511) produced in HEK293T cell.

Target HNRNPH1

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG1

Clonality Monoclonal


MW 49 kDa

Uniprot ID [P31943](#)

NCBI [NM_005520](#)

Dilution Range WB 1:500~2000, IHC 1:500

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...