


Product Datasheet

MASA (ENOPH1) antibody (orb1327578)

Description	Mouse monoclonal antibody to MASA (ENOPH1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human ENOPH1 (NP_067027) produced in E.coli.
Target	ENOPH1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	28.8 kDa
Uniprot ID	Q9UHY7
NCBI	NM_021204
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...