

Product Datasheet

CPLANE2 antibody (orb1327633)

Description	Mouse monoclonal antibody to CPLANE2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human RSG1 (NP_112169) produced in E.coli.
Target	CPLANE2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	28.3 kDa
Uniprot ID	Q9BU20
NCBI	NM_030907
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.

