


Product Datasheet

LINC02145 antibody (orb1327655)

Description	Mouse monoclonal antibody to LINC02145
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human FLJ33360 (NP_001001702) produced in E.coli.
Target	FLJ33360
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	14.2 kDa
NCBI	NM_001001702
Dilution Range	WB 1:500, IHC 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...