

Product Datasheet

NAT1 antibody (orb1327688)

Description	Mouse monoclonal antibody to NAT1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human NAT1 (NP_000653) produced in E.coli.
Target	NAT1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P18440
NCBI	NM_000662
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.

