


Product Datasheet

BAG5 antibody (orb1327701)

Description	Mouse monoclonal antibody to BAG5
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human BAG5 (NP_004864) produced in 293T.
Target	BAG5
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	51 kDa
Uniprot ID	Q9UL15
NCBI	NM_004873
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...