

Product Datasheet

SERPINB8 antibody (orb1327781)

Description	Mouse monoclonal antibody to SERPINB8
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 169-374 of human SERPINB8 (NP_001027018) produced in E.coli.
Target	SERPINB8
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	27.5 kDa
Uniprot ID	P50452
NCBI	NM_001031848
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.

