


Product Datasheet

CDY2B antibody (orb1327817)

Description	Mouse monoclonal antibody to CDY2B
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human CDY2B (NP_001001722) produced in E.coli.
Target	CDY2B
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	60.3 kDa
Uniprot ID	Q9Y6F7
NCBI	NM_001001722
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...