


Product Datasheet

DHRS4L2 antibody (orb1327945)

Description	Mouse monoclonal antibody to DHRS4L2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human DHRS4L2 (NP_932349) produced in E.coli.
Target	DHRS4L2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	Q6PKH6
NCBI	NM_198083
Dilution Range	WB 1:500~2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...