


Product Datasheet

PTCD3 antibody (orb1327954)

Description	Mouse monoclonal antibody to PTCD3
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-275 of human PTCD3 (NP_060422) produced in E.coli.
Target	PTCD3
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	78.4 kDa
Uniprot ID	Q96EY7
NCBI	NM_017952
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...